



Obstetrical Tools of the Trade

Guenther's Embryotomy Knife



Thygesen's Fetatome



Utrecht Fetatome



Cleaning Brush for fetatome



OB Saw Wire (Gigli Wire)



Lyss OB Handles



Sutter Style
OB Handle



OB Wire
Handle

OB Wire Guide (“wire passer”)



OB Wire Drop Guide




Krey Hook



Vienna Double Blunt Eye Hook



Moore's OB Chain Handle, OB Chain, and proper placement of chains on fetlock



First loop is placed proximal to the fetlock, followed by a half hitch distal to the fetlock



Calf ("Head") Snare



Torsion Rod



A photograph of a Hercules Fetal Extractor in a clinical setting. The device is a long, silver, corrugated metal pole with a red handle and a chain attached to a horizontal bar above. A person's arm is visible holding the chain. The background shows a room with metal railings, a door, and a blue light fixture.

Hercules Fetal Extractor